


# soft**Flow** Homecare Applicators

Patient interface for homecare application in adult and adolescent patients





- ✓ High wearing comfort due to individual fit, simple handling and pleasant skin tolerance
- ✓ Soft and smooth thanks to medical silicone
- ✓ Patented, integrated heating up to the nose prevents condensation
- ✓ 7 different applicator types for optimal comfort at any flow rate
- ✓ Optional interface for tracheotomy patients with flow rates up to 50 l/min



# softFlow Homecare Applicators

## TECHNICAL DATA:

| Type:                                                                                                  | Applicator Homecare Small                                                           | Applicator Homecare Small-Plus                                                      | Applicator Homecare Standard                                                        | Applicator Homecare Standard-Plus                                                   | Applicator Homecare Medium                                                           | Applicator Homecare Large                                                             | Tracheal Applicator Homecare                                                               |
|--------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------|
| Article number:                                                                                        | 40630101                                                                            | 40630104                                                                            | 40630102                                                                            | 40630105                                                                            | 40630110                                                                             | 40630113                                                                              | 40630119                                                                                   |
| Hygiene:                                                                                               | Single patient multiple use                                                         | Single patient multiple use                                                         | Single patient multiple use                                                         | Single patient multiple use                                                         | Single patient multiple use                                                          | Single patient multiple use                                                           | Single patient multiple use                                                                |
| Maximum flow rate:                                                                                     | 20 l/min                                                                            | 45 l/min                                                                            | 25 l/min                                                                            | 35 l/min                                                                            | 60 l/min                                                                             | 60 l/min                                                                              | 50 l/min                                                                                   |
| Operating life:                                                                                        | 720 therapy hours                                                                   | 720 therapy hours                                                                   | 720 therapy hours                                                                   | 720 therapy hours                                                                   | 720 therapy hours                                                                    | 720 therapy hours                                                                     | 720 therapy hours                                                                          |
| Distance NC (A):<br>   | 12,5 mm                                                                             | 12,5 mm                                                                             | 15 mm                                                                               | 15 mm                                                                               | 15 mm                                                                                | 16 mm                                                                                 | -                                                                                          |
| Length NC (B):<br>     | 11 mm                                                                               | 11 mm                                                                               | 13 mm                                                                               | 13 mm                                                                               | 13 mm                                                                                | 13 mm                                                                                 | -                                                                                          |
| Diameter NC (D):<br> | 3,2 mm                                                                              | 3,2 mm                                                                              | 4,0 mm                                                                              | 4,0 mm                                                                              | 4,0 mm                                                                               | 5,0 mm                                                                                |  22 mm |
| System connection:                                                                                     |  |  |  |  |  |  |       |
| Total length (L):                                                                                      | 225 cm                                                                              | 225 cm                                                                              | 225 cm                                                                              | 225 cm                                                                              | 225 cm                                                                               | 225 cm                                                                                | 190 cm                                                                                     |
| Supply tube length:<br>Inner diameter:<br>Outside diameter:                                            | 180 cm<br>6 mm<br>9 mm                                                              | 180 cm<br>8 mm<br>12 mm                                                             | 180 cm<br>6 mm<br>9 mm                                                              | 180 cm<br>8 mm<br>12 mm                                                             | 180 cm<br>8 mm<br>12 mm                                                              | 180 cm<br>8 mm<br>12 mm                                                               | 180 cm<br>7 mm<br>13 mm                                                                    |
| Loop (S) length:<br>Inner diameter:<br>Outside diameter:                                               | 2x45 mm<br>3,6 mm<br>5 mm                                                           | 2x45 mm<br>5 mm<br>7 mm                                                             | 2x45 mm<br>3,6 mm<br>5 mm                                                           | 2x45 mm<br>3,6 mm<br>5 mm                                                           | 2x45 mm<br>5 mm<br>7 mm                                                              | 2x45 mm<br>5 mm<br>7 mm                                                               | -                                                                                          |
| Heating:                                                                                               | max. 24 VA                                                                          | max. 24 VA                                                                          | max. 24 VA                                                                          | max. 24 VA                                                                          | max. 24 VA                                                                           | max. 24 VA                                                                            | max. 24 VA                                                                                 |
| CE Certificate:                                                                                        | CE 0297                                                                             | CE 0297                                                                             | CE 0297                                                                             | CE 0297                                                                             | CE 0297                                                                              | CE 0297                                                                               | CE 0297                                                                                    |
| Disposal:                                                                                              | Domestic waste                                                                      | Domestic waste                                                                      | Domestic waste                                                                      | Domestic waste                                                                      | Domestic waste                                                                       | Domestic waste                                                                        | Domestic waste                                                                             |
| Manufacturer:                                                                                          | TNI medical AG<br>Nürnberger Str. 74a<br>97076 Würzburg<br>Germany                  | TNI medical AG<br>Nürnberger Str. 74a<br>97076 Würzburg<br>Germany                  | TNI medical AG<br>Nürnberger Str. 74a<br>97076 Würzburg<br>Germany                  | TNI medical AG<br>Nürnberger Str. 74a<br>97076 Würzburg<br>Germany                  | TNI medical AG<br>Nürnberger Str. 74a<br>97076 Würzburg<br>Germany                   | TNI medical AG<br>Nürnberger Str. 74a<br>97076 Würzburg<br>Germany                    | TNI medical AG<br>Nürnberger Str. 74a<br>97076 Würzburg<br>Germany                         |

